























## SynergEyes, Inc. Symbol Glossary

---

Symbol	Symbol Title	Description of Symbol	Standard Reference*	Standard Title
<b>Symbols Derived from Standards</b>				
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC, and 98/79/EC.	ISO 15223-1, Symbol 5.1.1	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Authorized representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, Symbol 5.1.2	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1, Symbol 5.1.4	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1, Symbol 5.1.5	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, Symbol 5.1.6	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Serial number	Indicates the manufacturer's serial number so that a specific medical device can be identified.	ISO 15223-1, Symbol 5.1.7	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.	ISO 15223-1, Symbol 5.2.4	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Sterilized using steam or dry heat	Indicates a medical device that has been sterilized using steam or dry heat.	ISO 15223-1, Symbol 5.2.5	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements

Symbol	Symbol Title	Description of Symbol	Standard Reference*	Standard Title
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened. This symbol may also mean “Do not use if the product sterile barrier system or its packaging is compromised.”	ISO 15223-1, Symbol 5.2.8	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.	ISO 15223-1, Symbol 5.3.2	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Keep dry	Indicates a medical device that needs to be protected from moisture.	ISO 15223-1, Symbol 5.3.4	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1, Symbol 5.3.7	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1, Symbol 5.4.2	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1, Symbol 5.4.3	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1, Symbol 5.4.4	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
<b>Symbols Not Derived from Standards</b>				
	Prescription Only	Requires prescription in the United States.  Caution: Federal (USA) law restricts this device to sale	21 CFR 801.15	Part 801 Labeling, Section 801.15 Medical devices; prominence of required label statements; use of symbols in labeling

Symbol	Symbol Title	Description of Symbol	Standard Reference*	Standard Title
		by or on the order of a licensed eyecare practitioner.	21 CFR 801.109	Part 801 Labeling, Section 801.109 Prescription devices
	CE marking of conformity	Signifies European technical conformity	93/42/EEC Articles 4,11,12,17, Annex II	the requirements for accreditation and market surveillance relating to the marketing of products;
	Trash box, litter bin, rubbish bin	To indicate a receptacle for the acceptance of trash or litter or rubbish being thrown away	ISO 7001	Graphical symbols – public information symbols
	UV Blocking	Contact lens contains UV radiation blocking properties.	-	-
	Blue tinted lens	-	-	-
	Right Eye	-	-	-
	Left Eye	-	-	-
<b>Abbreviations</b>				
DIA	Diameter			
BC	Base Curve			
VLT	Vault			
SC	Skirt Curve			
PWR (P)	Diopter Power			
ADD	Add Power (Multifocal)			
ZN	Zone (Multifocal)			
LIFT	Lift			
CYL	Cylinder			
TAN	Tangent			
CT	Center Thickness			
MAT	Material			
DIAG	Diagnostic			

\* References to ISO 15223-1:2016 and BS EN ISO 15223-1:2016.