

SynergEyes, Inc. Symbol Glossary

Document Number: 70201

Revision: D

Release Date: 2025-10-07

ECR: 2096

Symbol	Symbol Title	Description of Symbol	Standard Reference*	Standard Title
Symbols Derived from Standards				
	Manufacturer	Indicates the medical device manufacturer	ISO 15223-1, Symbol 5.1.1	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Date of Manufacture Country of Manufacture	Indicates the date when the medical device was manufactured. Indicates that country of manufacture of products.	ISO 15223-1 Symbol 5.1.3 Symbol 5.1.11	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Medical Device (MD)	Indicates the item is a medical device	ISO 15223-1 Symbol 5.7.7	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Authorized representative	Indicates the authorized representative in the identified country or jurisdiction (e.g., European Community/European Union).	ISO 15223-1, Symbol 5.1.2	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1, Symbol 5.1.4	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1, Symbol 5.1.5	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, Symbol 5.1.6	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Serial number	Indicates the manufacturer's serial number so that a specific	ISO 15223-1, Symbol 5.1.7	Medical devices – Symbols to be used with medical device labels, labelling and information

SynergEyes, Inc. Symbol Glossary

Document Number: 70201

Revision: D

Release Date: 2025-10-07

ECR: 2096

Symbol	Symbol Title	Description of Symbol	Standard Reference*	Standard Title
		medical device can be identified.		to be supplied – Part 1: General requirements
	Unique device identifier (UDI)	Indicates a carrier that contains unique device identifier information	ISO 15223-1 Symbol 5.7.10	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Sterilized using irradiation Single sterile barrier system with protective packaging outside	Indicates a medical device that has been sterilized using irradiation. Indicates a single sterile barrier system with protective packaging outside.	ISO 15223-1, Symbol 5.2.4 Symbol 5.2.14	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened. This symbol may also mean “Do not use if the product sterile barrier system or its packaging is compromised.”	ISO 15223-1, Symbol 5.2.8	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1, Symbol 5.3.7	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Single patient multiple use	Indicates a medical device that may be used multiple times (multiple procedures) on a single patient.	ISO 15223-1, Symbol 5.4.12	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1, Symbol 5.4.2	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	Consult instructions for use or consult	Indicates the need for the user to consult the instructions for use.	ISO 15223-1, Symbol 5.4.3	Medical devices – Symbols to be used with medical device labels, labelling and information

SynergEyes, Inc. Symbol Glossary

Document Number: 70201

Revision: D

Release Date: 2025-10-07

ECR: 2096

Symbol	Symbol Title	Description of Symbol	Standard Reference*	Standard Title
	electronic instructions for use			to be supplied – Part 1: General requirements
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1, Symbol 5.4.4	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
Symbols Not Derived from Standards				
 	Prescription Only	Requires prescription in the United States. Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed eyecare practitioner.	21 CFR 801.15	Part 801 Labeling, Section 801.15 Medical devices; prominence of required label statements; use of symbols in labeling
			21 CFR 801.109	Part 801 Labeling, Section 801.109 Prescription devices
	Diagnostic Lens	Symbol indicating that the lens is intended for fitting and is not a prescription lens.	-	-
	CE marking of conformity	Signifies European technical conformity.	93/42/EEC Articles 4,11,12,17, Annex II	The requirements for accreditation and market surveillance relating to the marketing of products;
	Do not use with water, Avoid contact with water	To indicate that the contact lenses and contact lens cases should not be exposed to water, and to indicate that contact lenses should not be worn in the water.	-	-
	UV Blocking	Contact lens contains UV radiation blocking properties.	-	-
	Blue tinted lens	-	-	-

SynergEyes, Inc. Symbol Glossary

Document Number: 70201

Revision: D

Release Date: 2025-10-07

ECR: 2096

Symbol	Symbol Title	Description of Symbol	Standard Reference*	Standard Title
	Green tinted lens			
	Violet tinted lens			
	Right Eye	-	-	-
	Left Eye	-	-	-
Abbreviations				
DIA	Diameter			
BC	Base Curve			
VLT	Vault			
SC	Skirt Curve			
PWR (P)	Diopter Power			
ADD	Add Power (Multifocal)			
ZN	Zone (Multifocal)			
LIFT	Lift			
CYL	Cylinder			
TAN	Tangent			
CT	Center Thickness			
MAT	Material			
DIAG	Diagnostic			
AC	Aberration Control			
LC	Limbal Clearance Factor			

* References to ISO 15223-1:2021+A1:2025 and BS EN ISO 15223-1:2021+A1:2025.